



Verza™

A New Frontier in Ultrasound Guidance

Designed with extensive physician input, the Verza™ Guidance System delivers unmatched clinical versatility by being the first commercial guide to provide five angles as well as the industry's largest range of interventional devices (14 gauges) for precise instrument placement.

Verza's™ unique design includes VerzaLink™, a direct attach locating feature built into the transducer, offering simplified attachment to the transducer for optimal performance.

Features:

- Reduced number of blind zones and unguided needle length leads to confident outcomes.¹
- Increased channel length with an improved accuracy of 25%.¹
- Reduced technique variability provides a shorter learning curve and reduced procedure time.
- Improved user experience due to bracket-free design (direct transducer attachment), which enhances transducer handling.

Quick-release function, allowing easy detachment of needle.

Five-angle approach for improved anatomical access and an increased gauge range.



¹ CIVCO reporting data on file. All comparisons based on CIVCO with conventional brackets versus Verza™ direct attachment guidance system



Ordering Information

Verza™ Needle Guide Disposable Replacement Kits

Ref no.	Descriptions	Qty/box
610-1500-24	Sterile Verza™ with 14 x 91.5cm telescopically-folded CIV-Flex™ cover	24
610-1501-24	Sterile Verza™ with 14 x 147cm telescopically-folded CIV-Flex™ cover	24
610-1502-24	Sterile Verza™ with 8.9 x 91.5cm telescopically-folded CIV-Flex™ cover	24
610-1504-24	Sterile Verza™ with 14 x 147cm telescopically-folded CIV-Flex™ cover with PEELSafe Advantage	24

MM 02-1000-03 ENG Rev. 02.

Contact Details, Nordics

customer@mermaidmedical.com
mermaidmedical.com

Contact Details, NL

info@mermaidmedical.nl
mermaidmedical.com

Contact Details, UK

customer@mermaidmedical.co.uk ·
mermaidmedical.com